

# Pump In Style® with maxFLOW™ Double Electric Breast Pump

Now Available at Pharmacare

Pump in Style is a double electric breast pump developed for moms who pump several times a day. Trusted performance from the #1 recommended breast pump brand<sup>1</sup>, Pump in Style is designed to make breastfeeding easier, so you can feed your baby with amazing breast milk for longer. Featuring MaxFlow™ technology and 2-Phase Expression® means more milk in less time.



## Benefits



**Hospital performance**



Few parts – **easy to clean and assemble**



Removes **11.8% more milk** per minute<sup>2</sup>



**2-Phase Expression Technology** mimics baby's natural sucking rhythm



**Closed system** at the kit prevents breast milk from entering the tubing\*



**PersonalFit™ PLUS** breast shields with comfortable soft rim and oval shape for a better fit<sup>2</sup>



Intuitive controls for **ease of use**

**Efficient:** Advanced motor generates vacuum with MaxFlow micro-vibrations for effective breast milk expression

**Designed For Daily Use:** Double electric breast pump designed for moms who pump several times a day

**Safe For You and Baby:** All parts that touch breast milk are made without BPA

**Convenience:** Includes a battery pack for convenient pumping away from the power outlet\*\*

<sup>1</sup>IQVIA ProVoice Survey; September 2018 – August 2019.

<sup>2</sup>Clinical studies, PersonalFit Flex™ compared to PersonalFit™ data on file (NCT02496429, NCT02492139, NCT02719548, NCT03091985).

\*Using other connectors can cause irreversible damage to your breast pump.

\*\*Batteries not included

## Included Components:



## Product Specifications:

Item #	Description	Qty Per Case	Case Dimensions (HxWxD) Inches	Case Weight lbs	Qty/Case (each)
101041360	NEW Pump in Style Breast Pump	1	4.65 x 11.73 x 9.72	3.28	1

## Technical Specifications:

Vacuum (mmHg)	-50 to -295mmHg Maximum vacuum - expression (double and single pumping).
CPM	54 to 120

**Important:** Plastic bottles and component parts become brittle when frozen and may break when dropped. Also, bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened, or knocked over. Take appropriate care in handling bottles and components. Do not use the breast milk if bottles or components become damaged.

All components that come in contact with mother's milk are made without BPA and meet US Consumer Product Safety Improvement Act (CPSIA 2008), CA AB 1108 requirements for DEHP and other phthalates for all components that contact mother's milk.

Warranted to the original participant to be free from defects in material and workmanship for the period of one year for pump mechanism (90 days for parts and accessories).

[www.pharmacarehawaii.com/breast-pumps](http://www.pharmacarehawaii.com/breast-pumps)

**To order or for more information, contact: Pharmacare's Breast Pump Hotline 808-840-4199.**